

## Gynaec Examination cum Delivery table, Mechanical

Product Model : **HF302**



**Category: Obstetric Labour Tables, hospital furniture, General Medical Devices,**

### Specification :

- MS Tubular frame, pre-treated and epoxy powder coated.
- Fixed height
- 3 section cushioned top, framework mounted on rubber shoes.
- Padded knee crutches adjustable.
- Knockdown Construction, Sides of this section are fit with vertical handgrips.
- Optional:- Lithotomy Pads / Lithotomy Striups.
- Supplied with S.S Basin
- Available in Stainless Steel (Code):- **HF302.SS**

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com