

7.0mm Self Drilling Cannulated Cancellous Screw, Hexagonal Socket-32mm Threaded

Product Model : **2207**



Category: Cannulated Cancellous Screw, Orthopaedic Implants,


Specification :

- 7.0mm Self Drilling Cannulated Cancellous Screw, Hexagonal Socket-32mm Threaded

Item Code (Titanium)	Length(mm)
Ti-2207.50	50 mm
Ti-2207.55	55 mm
Ti-2207.60	60 mm
Ti-2207.65	65 mm
Ti-2207.70	70 mm
Ti-2207.75	75 mm
Ti-2207.80	80 mm
Ti-2207.85	85 mm

	Ti-2207.90	90 mm
	Ti-2207.95	95 mm
	Ti-2207.100	100 mm
	Ti-2207.105	105 mm
	Ti-2207.110	110 mm
	Ti-2207.115	115 mm

Legacy of Excellence in Health Care..

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com