

6.5mm Cannulated Cancellous Screw, Hexagonal Socket-Fully Threaded

Product Model : **2211**



Category: Cannulated Cancellous Screw, Orthopaedic Implants,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Available in S.S & Titanium Material.

Thread Dia	6.5 mm
Head Dia	8.0 mm
Core Dia	4.5 mm
Lengths	50-130 mm
Guide wire used	1.8 mm
Cannulated Drill Bit	4.5 mm
Thread Length	Fully threaded