

## 6.5mm Cannulated Cancellous Screw, Hexagonal Socket-Fully Threaded

Product Model : **2211**




**Category: Cannulated Cancellous Screw, Orthopaedic Implants,**

**Specification :**

- Available in S.S & Titanium Material.

<b>Thread Dia</b>	6.5 mm
<b>Head Dia</b>	8.0 mm
<b>Core Dia</b>	4.5 mm
<b>Lengths</b>	50-130 mm
<b>Guide wire used</b>	1.8 mm
<b>Cannulated Drill Bit</b>	4.5 mm

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)