

3.5mm Small Cannulated Cancellous Screw, Partially Threaded

Product Model : **2201**



Category: Cannulated Cancellous Screw, Orthopaedic Implants,

Legacy of Excellence in Healthcare...

Description

- Available in S.S & Titanium Material.

Technical Specification	
Thread Dia	3.5mm
Drill Bit for Threaded Hole	2.5mm
Core Dia.	1.9mm
Hexagonal Socket Width	2.5mm
Head Dia.	6.0mm
Thread Length	Partially
Lengths	10mm to 60mm