

## Carbon Steel Surgical Blade

Product Model : **SB01**



**Category: Surgical Blades, Medical Equipment,**

### Specification :

- Surgical blades are available in the following sizes 10, 10A, 11, 11R, 12, 15, 15D, 18, 20, 21, 22, 22A, 23, 24, 25 & 36.
- They are made of either Carbon steel or Stainless steel (sterile or non-sterile) singly packed in peel open pouches with a VCI protective inner liner.
- A convenient number of individual pouches are boxed in a unit carton.
- The surgical blades are sterilized by GAMMA RADIATION

Carbon Steel Surgical Blade							
Size	10	10A	11	11P	12	15	16
Item Code	C2	C3	C4	C5	C6	C9	C12
Blade Handle No.	3	3	3	3	3	3	3

  

Size	18	19	20	21	22	23	24
------	----	----	----	----	----	----	----

Item Code	C14	C15	C16	C17	C18	C20	C21
Blade Handle No.	4	4	4	4	4	4	4

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com