

## Carbon Steel Surgical Blade

Product Model : **SB01**



**Category: Surgical Blades, General Medical Devices,**

### Specification :

- Surgical blades are available in the following sizes 10, 10A, 11, 11R, 12, 15, 15D, 18, 20, 21, 22, 22A, 23, 24, 25 & 36.
- They are made of either Carbon steel or Stainless steel (sterile or non-sterile) singly packed in peel open pouches with a VCI protective inner liner.
- A convenient number of individual pouches are boxed in a unit carton.
- The surgical blades are sterilized by GAMMA RADIATION

Carbon Steel Surgical Blade							
<b>Size</b>	10	10A	11	11P	12	15	16
<b>Item Code</b>	C2	C3	C4	C5	C6	C9	C12
<b>Blade Handle No.</b>	3	3	3	3	3	3	3
<b>Size</b>	18	19	20	21	22	23	24

Item Code	C14	C15	C16	C17	C18	C20	C21
Blade Handle No.	4	4	4	4	4	4	4

*Legacy of Excellence in Health Care..*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com