

Disposable Isolation Gown, SMS

Product Model : **Default Model No**




Front View

Category: Non Woven, Medical Disposables, Medical Equipment, Covid-19,

Summary :

Discover high-quality Disposable Isolation Gowns (SMS) designed to protect healthcare professionals and patients from cross-contamination in hospitals, clinics, and medical facilities. Made from durable, breathable, and fluid-resistant SMS material, these gowns provide comfortable fit, reliable protection, and easy disposal. Ideal for infection control, medical procedures, and daily patient care, ensuring hygiene and safety in every setting.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- SMS non-woven or SMMS non-woven fabric
- GSM Available: 25 to 70g
- Size: L, XL, XXL, XXXL, or as Your Request
- Breathable, dust-proof, waterproof, filter bacteria high density.
- Adjustable tie type waistline.
- Special neckline heightening.
- Knitted Cotton cuffs.

Code	Sizes	Size Details (in cms)				
		Length	Chest	Sleeve Length	Sleeve Opening	Ankle Opening
DISPC22.1	Small	165	120	84	18	24
DISPC22.2	Medium	169	125	86	18	24
DISPC22.3	Large	173	130	90	18	24
DISPC22.4	X-Large	178	135	93	18	24
DISPC22.5	XX-Large	181	140	96	18	24
DISPC22.6	XXX-Large	188	145	99	18	24
		Error:- ±2	Error:- ±2	Error:- ±2	Error:- ±2	Error:- ±2