

## Shoe Cover - PE (Polyethylene), Disposable?

Product Model : **DISP06**



**Category: Covid-19, Medical Equipment, Medical Disposables, Non Woven,**

### Summary :

Shoe Cover, Polyethylene Single Use, Disposable. These are single use shoe covers that can be worn and removed easily. Our disposable plastic shoe cover is used for a variety of applications in medical industries like hospitals and pharmaceuticals. We provide a barrier against possible exposure to airborne organisms or contact with a contaminated environment.

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Disposable Shoe Cover
- Box. of 100pcs. (50 pairs)
- Made from polyethylene (LDPE) film.
- Elastic bands around the shoe give secure comfort fit.
- Good resistance to liquids
- Gaze:- 90
- Size:- 30cms x 16.5cms
- Available in Green & Blue Color.
- Master Cartoon Dimension (inches):- 22? x 17? x 19?
- Qty per Master Cartoon:- 5000 pcs. / 50 Box of 100pcs.