

## Shoe Cover - PE (Polyethylene), Disposable?

Product Model : **DISP06**



**Category: Covid-19, General Medical Devices, Medical Disposables, Non Woven,**

### Summary :

Shoe Cover, Polyethylene Single Use, Disposable. These are single use shoe covers that can be worn and removed easily. Our disposable plastic shoe cover is used for a variety of applications in medical industries like hospitals and pharmaceuticals. We provide a barrier against possible exposure to airborne organisms or contact with a contaminated environment.

### Specification :

- Disposable Shoe Cover
- Box. of 100pcs. (50 pairs)
- Made from polyethylene (LDPE) film.
- Elastic bands around the shoe give secure comfort fit.
- Good resistance to liquids
- Gaze:- 90
- Size:- 30cms x 16.5cms
- Available in Green & Blue Color.
- Master Cartoon Dimension (inches):- 22? x 17? x 19?
- Qty per Master Cartoon:- 5000 pcs. / 50 Box of 100pcs.

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com