

## Anaesthesia Machine Major "A"

Product Model : **AN-161**



**Category: Anaesthesia Products & Equipment, Anaesthesia Machine,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Capacious drawer for keeping anaesthesia kit/accessories.
- Extended platform at rear for two 10ltrs. (Water Cap.) Cylinders.
- Large roller bearings, antistatic castors, for easy mobility and brakes.
- Top tray provided at eye level for monitoring equipments.
- Hook for breathing circuit.
- Clean stainless steel top.
- Provision to fit ventilator.
- Emergency drug tray fitted at convenient place Unitised rota meter block for oxygen & nitrous oxide with fine adjustment controls & colour coated knob and luminous back plate. Oxygen flow meter calibrated from 100cc/min to 10 ltrs./ Minute & Nitrous oxide flow meter from 200cc/mm to 12 liters/min.
- One gold mantype vaporizer for halothene, provision to mount three vaporizer at a time.
- Non return cum blow off valve in circuit ensures that pressure does not build-up beyond 200cm. Water column & minimises risk of back flow of gases.
- 60 psl outlet for driving ventilator.
- Change over system open circuit or close circuit.
- Emergency Oxygen flush valve at table level to give 35-60 ltrs/min for emergencies.
- Cylinder yokes: pin indexed yokes two each of oxygen & nitrous oxide.
- 63 mm diameter, colour coded gauges fitted at convenient angle for better visibility & clarity, one gauge for O<sub>2</sub> & one gauge for
- N<sub>2</sub>O shown the line pressure
- High efficiency regulators for reliability, two each for O<sub>2</sub> & N<sub>2</sub>O
- One Maggills or Brain Circuit
- One B.P Apparatus.