

## Anaesthesia Machine Mini "T"

Product Model : **AN-163**



**Category: Anaesthesia Products & Equipment, Anaesthesia Machine,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Construction: Trolley made of tubular rigid steel section electrostatically powder coated or stainless steel
- Suitable drawer for keeping anaesthesia kit/accessories
- Extended platform at rear for one 10 ltrs. (Water capacity) cylinder
- Small roller bearings, antistatic castors, for easy mobility and brakes
- Top tray provided at eye level for monitoring equipments
- Hook for breathing circuit
- Wipe clean stainless steel top
- Unitized rotameter block for oxygen & nitrous oxide with fine adjustment controls & colour coated knob and luminous back plate. Oxygen flow meter calibrated from 100cc/min to 10 ltrs. Minute & Nitrous oxide flow meter from 200/min to 12 liters/min down Stream receipt of oxygen
- One gold man type vaporizer for halothene
- Non return cum blow off valve in circuit ensures that pressure does not build-up beyond 200cm. water column & minimises risk of back flow of gases.
- 60 psl outlet for driving ventilator
- Emergency Oxygen flush valve at table level to give 35-60 ltrs./min for emergencies
- Cylinder yokes: Pin indexed yokes, two each of oxygen & nitrous oxide
- 50mm diameter, colour coded gauges fitted at convenient angle for better visibility & clarity
- High efficiency regulators for reliability, one each for O<sub>2</sub> & N<sub>2</sub>o
- One Maggill's or Brain circuit