

Green Max “SPR” Blades - Fiber Optic (4.6mm Light Guide)

Product Model : **LSC21**



Category: Medical Equipment, Laryngoscope Set & Spares,

Summary :

GWS Fiber Optic Laryngoscopes deliver superior illumination for safer intubations. **Reusable Fiber Optic blades** are manufactured with 100% high quality Stainless Steel to ensure durability and long life. Each blade is MATT finished to reduce glare and enhance ease of intubation for optimal patient safety.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Blades with non-replaceable fibre optics.
- All blades are made of premium grade stainless steel in matt finish.
- Fibre optics Ø 4.6mm / 6800 single fibres minimum.
- Environment friendly, economical & cost effective.
- Can be treated in autoclave up to 134° - approx. 4000 times.
- In accordance with DIN ISO 7376

Size	Overall Length	Blade Length	Distal Width	Item Code
MAC0	75 mm	58 mm	9 mm	LSC21.1
MAC1	92 mm	72 mm	9 mm	LSC21.2
MAC2	101 mm	93 mm	13 mm	LSC21.3
MAC3	133 mm	110 mm	14.5 mm	LSC21.4
MAC4	155 mm	130 mm	15.5 mm	LSC21.5
MAC5	172 mm	150 mm	15.5 mm	LSC21.6
MILL00	65 mm	44 mm	10 mm	LSC21.7
MILL0	75 mm	54 mm	10 mm	LSC21.8
MILL1	102 mm	79 mm	10 mm	LSC21.9
MILL2	155 mm	131 mm	13 mm	LSC21.10
MILL3	195 mm	172 mm	13 mm	LSC21.11
MILL4	205 mm	182 mm	17 mm	LSC21.12