


Flexitip Optic LED Balde

Product Model : **LSC13**



Category: Medical Equipment, Laryngoscope Set & Spares,

Summary :

The Forte **GWS Flitip blade** is specially designed to improve visualization of the larynx during intubations. With a blade tip adjustable through 80degrees the Forte Flitip blade can simplify elevation of the epiglottis for many conditions, including reduced thyromental distance, short or immobilized neck, obesity and swollen tissue of the larynx.

Forte Flitip laryngoscope blades are ideal for EMS, Emergency Room, and any other setting with unpredictable patient airway conditions.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- LED technology provides up to 40% more light than conventional halogen bulbs.
- MATT finished blade surface reduces glare to enhance ease of intubation.
- Compatible with all ISO 7376-1 Conventional handles.

Size	Item Code
MAC2	LSC13.1
MAC3	LSC13.2
MAC3.5	LSC13.3
MAC4	LSC13.4
MAC4.5	LSC13.5