

Shallow Tray, Stainless Steel

Product Model : **HH-18**



Category: General Medical Devices, Stainless Steel Hospital Hollowware,

Specification :

Shallow Tray (Stainless Steel)			
Approx. Size		Quality / Code	
In Inches	In mm	Deluxe Quality/ Thickness: 0.5mm	Super Deluxe Quality / Thickness: 0.5mm/ 304 Grade
8" x 6"	200 x 150 x 25	HH-18.1.1	HH-18.2.1
9" x 6"	230 x 155 x 25	HH-18.1.2	HH-18.2.2
10" x 8"	255 x 200 x 25	HH-18.1.3	HH-18.2.3

11" x 7"	280 x 180 x 25	HH-18.1.4	HH-18.2.4
12" x 8"	300 x 205 x 25	HH-18.1.5	HH-18.2.5
12" x 10"	305 x 255 x 25	HH-18.1.6	HH-18.2.6
14" x 10"	355 x 255 x 25	HH-18.1.7	HH-18.2.7
15" x 12"	380 x 305 x 25	HH-18.1.8	HH-18.2.8
18" x 12"	455 x 305 x 25	HH-18.1.9	HH-18.2.9

Legacy of Excellence in Health Care..



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com