

## Contoured L.S Belt - Primo

Product Model : **RHA514**



**Category: General Medical Devices, Rehabilitation Supports, Back Support,**

### Specification :


**Features:** Highly porous elastic webbing improves ventilation and comfort. Anatomically contoured back splints take the shape of the body and provide comfortable immobilization and help maintain correct posture. The Elastic Contoured LS Support provides ideal back support for back ache. Heat resistant rubber with the high modulus of elasticity improves the life of the belt. Double pull mechanism. Well cushioned back pad. Anatomical pre-shaped splints.

### Application:

- Degenerative Disc Disease.
- Post Dichotomy Care
- Sacroiliac mechanical due to poor posture & effectiveness for lower back ache patients.
- Dislocation of vertebrae.
- Low back pain
- Post back injury and fracture.
- Comfort during long working hours and travel.

Contoured L.S Belt - Primo		
Size	Size Inches	Item Code
Small	28-32	RHA514.1
Medium	34-36	RHA514.2
Large	38-40	RHA514.3
X-Large	42-44	RHA514.4
XX-Large	46-48	RHA514.5

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)