

## **L.S (Lumbar Sacro) Belt, Primo**

Product Model : **RHA512**



**Category: General Medical Devices, Rehabilitation Supports, Back Support,**

### **Specification :**

**Features:** Elastic supports on both sides for snug fit. Front and back panels are made of soft foam for comfort. Hook and loop closures for easy application and removal. Two rigid aluminium bars and two plastic bars for reinforcement. Double stitching on the edges for durability


### **Application:**

- Chronic low back pain & strain.
- Osteoporotic pain in the lumbosacral region.
- Degenerative Lumbar Disc Disease.
- Post surgical lumbo sacral immobilization.

Sacro Lumbar (L.S) Belt Primo	
Sizes	Item Code

Small	RHA512.1
Medium	RHA512.2
Large	RHA512.3
X Large	RHA512.4
XX Large	RHA512.5

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)