

Sacro Lumbar Brace (Belt)

Product Model : **RHA513**



Category: Medical Equipment, Rehabilitation Supports, Back Support,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- **Features:**

Made of tempered steel strips, it immobilizes the lumbar and sacral regions effectively. Fully elasticated to add extra pressure and support to the sacral and lumbar regions.

- **Applications:**

Degenerative lumbar disc disease. Osteoporotic pain of the lumbar spine. Post discectomy symptoms of the lumbar spine and lumbosacral strain.

- To relax sacrospinalis muscles to modify pelvic inclination and lumbar spine and make the patient constantly aware of the back pain.
-

Sacro Lumbar Brace (Belt)			
Size Inches	Quality / Codes		
	Super Fine Double Lock	Super Fine	Fine
28"	RHA513.1.1	RHA513.2.1	RHA513.3.1
30"	RHA513.1.2	RHA513.2.2	RHA513.3.2
32"	RHA513.1.3	RHA513.2.3	RHA513.3.3
34"	RHA513.1.4	RHA513.2.4	RHA513.3.4
36"	RHA513.1.5	RHA513.2.5	RHA513.3.5
38"	RHA513.1.6	RHA513.2.6	RHA513.3.6
40"	RHA513.1.7	RHA513.2.7	RHA513.3.7
42"	RHA513.1.8	RHA513.2.8	RHA513.3.8
44"	RHA513.1.9	RHA513.2.9	RHA513.3.9