

## Taylor Brace

Product Model : **RHA516**



**Category: Medical Equipment, Rehabilitation Supports, Back Support,**

### Summary :

Introducing the Taylor brace – a long spinal brace that provides excellent orthopedic support. This belt is ideal for bulk manufacturers and exporters globally, supplying top-quality orthopedic support to satisfied customers worldwide. The GWS brand guarantees excellence in products and service, supplying you with a reliable supplier of high-quality spinal braces.

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

**Features:**

Vertical as well as horizontal malleable aluminium stays to provide immobilization. Light weight splints with ethafoam cushioning improves comfort and patient compliance. Clavicle straps correct the shoulders posture while ensuring better fit. Double elastic panel allows easy application of the brace.

**Application:**

- Cox of spine.
- Osteoporosis.
- Mild fractures of the spine.

<b>Taylor Brace</b>		
<b>Size</b>	<b>Quality / Code</b>	
	<b>Regular</b>	<b>Deluxe</b>
Small	RHA516.1	RHA517.1
Medium	RHA516.2	RHA517.2
Large	RHA516.3	RHA517.3
X-Large	RHA516.4	RHA517.4