

## Taylor Brace

Product Model : **RHA516**



**Category: General Medical Devices, Rehabilitation Supports, Back Support,**

### Summary :

Introducing the Taylor brace – a long spinal brace that provides excellent orthopedic support. This belt is ideal for bulk manufacturers and exporters globally, supplying top-quality orthopedic supports to satisfied customers worldwide. The GWS brand guarantees excellence in products and service, supplying you with a reliable supplier of high-quality spinal braces.

### Specification :

**Features:** Vertical as well as horizontal malleable aluminium stays to provide immobilization. Light weight splints with ethafoam cushioning improves comfort and patient compliance. Clavicle straps correct the shoulders posture while ensuring better fit. Double elastic panel allows easy application of the brace.

### Application:

- Cox of spine.
- Osteoporosis.
- Mild fractures of the spine.

<b>Taylor Brace</b>		
<b>Size</b>	<b>Quality / Code</b>	
	<b>Regular</b>	<b>Deluxe</b>
Small	RHA516.1	RHA517.1
Medium	RHA516.2	RHA517.2
Large	RHA516.3	RHA517.3
X-Large	RHA516.4	RHA517.4

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com