

## Ambulance Collar

Product Model : **RHA529**



**Category: General Medical Devices, Rehabilitation Supports, Cervical Support,**

### Summary :


The Ambulance Collar is an essential piece of equipment in the event of a neck injury. This cervical collar supports and stabilizes the neck, preventing further damage and helping to ensure a speedy recovery. The Ambulance Collar is made from high-quality materials and is available at a competitive price.

### Specification :

- **Features:** Single Piece Flat Design. Semi Rigid Plastic, Padded With Sturdy Foam Provides Comfortable Immobilization. It Is Used For Temporary Cervical Immobilization For Accidental Victims During The Transportation & Diagnostic Studies.
- **Application:** It Is Used As Temporary Support To Neck In Accidental Conditions. It Provides Rigid Support To Limit The Neck Movement Or To Prevent Harmful Neck Movements.

<b>Ambulance Collar</b>	
Size	Item Code
Small	RHA529.1
Medium	RHA529.2
Large	RHA529.3

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmmed.com](mailto:info@gwsmmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmmed.com](http://www.gwsmmed.com)