

## Cervical Collar Hard Adjustable

Product Model : **RHA528**




**Category: Medical Equipment, Rehabilitation Supports, Cervical Support,**

### Summary :

Looking for a sturdy cervical collar that can be easily adjusted? Look no further than the Cervical Collar Hard Adjustable from GWS Surgicals LLP. This hard collar is made of high-quality materials and can be easily adjusted to provide the perfect level of support and comfort. Whether you're looking for an orthopedic support for daily use or a bulk manufacturer for export, GWS Surgicals LLP has you covered.

### Legacy of Excellence in Healthcare...

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- **Features:**

Made of light high density polythene sheet. Proper foam padding for comfortable wear. Double layer for height adjustment. Holes for ventilation. Velcro closure for easy application

- **Applications:**

Cervical spine immobilization in pre & post surgical cases.

Cervical spondylosis.

Cervical Collar Hard Adjustable	
Size	Item Code
Small	RHA528.1
Medium	RHA528.2
Large	RHA528.3
X-Large	RHA528.4