


Cervical Traction Kit

Product Model : **RHA533**




Category: Medical Equipment, Rehabilitation Supports, Cervical Support,


Summary :

Looking for a reliable cervical traction kit? Look no further than GWS Surgicals LLP. Our cervical traction kit is made from the highest quality materials and is backed by a 100% satisfaction guarantee. With our traction kit, you can enjoy relief from neck pain and other orthopedic issues. Send your bulk query today!

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com

 **Corp Office:** 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

Features:

Padding on the head halter to give comfortable wear. The spreader bar and pulley system enabled it to be attached to the bed. Heavy-duty weight bag with separate graduations for water & sand. Heavy-duty string is provided in the kit.

Applications:

Cervical spine traction traction in the sleeping position. Cervical spine dislocation. For continuous traction upon medical practitioner recommendation.

Size:

Universal

Cervical Traction Kit	
Item Code	Item / Models
RHA533.1	Cervical Traction Head Holder
RHA533.2	Cervical Traction Spreader Bar
RHA533.3	Traction Sand Bag
RHA533.4	Cervical Traction Kit (Sleeping Complete Set)