

Knee Brace Hinged

Product Model : **RHA503**




Category: Medical Equipment, Rehabilitation Supports, Knee/ Calf & Ankle Support's,


Summary :

Looking for a quality knee brace that will provide you with the support you need? Look no further than GWS Brand Knee Brace Hinged. This open patella brace is made of neoprene and is hinged to provide extra support and stability. Plus, it's available at a great price so you can get the relief you need without breaking the bank.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com

 **Corp Office:** 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

Hinged Knee Brace with Patella Stabilizer and knee protection pad for patella displacement prevention

Features:

Open patella with felt buttress pad relieves pressure on patella. Multidirectional neoprene circle pad around the patella provides padded protection and stabilization to prevent patella displacement. Two medial and lateral metal stays provide additional support and stability for the knee. Two hook and loop straps ensure comfortable fit and stabilization during activities. Helps relieve pain from tendonitis and arthritis. Retains body heat to the warm knee joint and increases blood circulation. Fit either the left or right Knee. Material: Neoprene, Velcro, Metal, Polyamide.

Application:

- Post ACL Ligament Tear Grade 1
- Post MCL/ LCL Ligament Tear Grade 1
- Post Operation Support
- Post Meniscal/ Cartilage Damage Grade 1

Knee Brace Hinged - Neo	
Sizes	Item Code
Small	RHA503.1
Medium	RHA503.2
Large	RHA503.3
X-Large	RHA503.4