

## **Knee Immobilizer 14"/19"/22"**

Product Model : **RHA504 | RHA505 | RHA506**



**Category: Medical Equipment, Rehabilitation Supports, Knee/ Calf & Ankle Support's,**

### **Summary :**

GWS Brand! Knee Immobilizer is perfect for immobilizing the knee joint and providing support. It is made of high-quality materials and can be used for both long and short term immobilization. Our Knee Immobilizer is also available in a variety of sizes like 14" / 19" / 22" to ensure a perfect fit.

### **Specification :**

**Features:** Its Contoured posterior and sturdy bilateral metal stays, immobilizes the leg at 0°-20°. It also has a removable medial, lateral and posterior stays. It has comfortable washable, skin friendly and laminated cloth material. Its patellar opening releases patellar pressure ensuring comfort. Its five straps of hook & loop fasteners help in secure fitting and better immobilization.

### **Application:**

- Post Total Knee Replacement.
- Post Ligament Reconstruction.

**Color Option:** Black/ Grey/ Beige

<b>Knee Immobilizer / Brace 14"</b>			
Item Code	Size	Inches	Cms
RHA504.1	Small	13-14	33-36
RHA504.2	Medium	14-15	36-38
RHA504.3	Large	15-16	38-41
RHA504.4	X-Large	16-18	41-46
RHA504.5	XX-Large	18-20	46-51

<b>Knee Immobilizer / Brace 19"</b>			
Item Code	Size	Inches	Cms
RHA505.1	Small	13-14	33-36
RHA505.2	Medium	14-15	36-38
RHA505.3	Large	15-16	38-41
RHA505.4	X-Large	16-18	41-46
RHA505.5	XX-Large	18-20	46-51

<b>Knee Immobilizer / Brace 22"</b>			
Item Code	Size	Inches	Cms
RHA506.1	Small	13-14	33-36
RHA506.2	Medium	14-15	36-38
RHA506.3	Large	15-16	38-41
RHA506.4	X-Large	16-18	41-46
RHA506.5	XX-Large	18-20	46-51

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com