

Knee Immobilizer 14"/19"/22"

Product Model : **RHA504 | RHA505 | RHA506**




Category: Medical Equipment, Rehabilitation Supports, Knee/ Calf & Ankle Support's,

Summary :

Discover high-quality knee immobilizers designed to provide support, stability, and comfort for patients with knee injuries or post-surgery recovery. Made from durable materials with adjustable straps, these knee braces ensure proper alignment and restrict movement to aid healing. Ideal for hospitals, clinics, and home care, knee immobilizers promote safety, recovery, and effective patient care.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

Features:

Its Contoured posterior and sturdy bilateral metal stays, immobilizes the leg at 0°-20°. It also has a removable medial, lateral and posterior stays. It has comfortable washable, skin friendly and laminated cloth material. Its patellar opening releases patellar pressure ensuring comfort. Its five straps of hook & loop fasteners help in secure fitting and better immobilization.

Application:

- Post Total Knee Replacement.
- Post Ligament Reconstruction.
- Emergency Immobilization for Knee

Color Option:

Black/ Grey/ Beige

Knee Immobilizer / Brace 14"			
Item Code	Size	Inches	Cms
RHA504.1	Small	13-14	33-36
RHA504.2	Medium	14-15	36-38
RHA504.3	Large	15-16	38-41
RHA504.4	X-Large	16-18	41-46
RHA504.5	XX-Large	18-20	46-51

Knee Immobilizer / Brace 19"			
Item Code	Size	Inches	Cms
RHA505.1	Small	13-14	33-36
RHA505.2	Medium	14-15	36-38
RHA505.3	Large	15-16	38-41
RHA505.4	X-Large	16-18	41-46

RHA505.5	XX-Large	18-20	46-51
----------	----------	-------	-------

Knee Immobilizer / Brace 22”			
Item Code	Size	Inches	Cms
RHA506.1	Small	13-14	33-36
RHA506.2	Medium	14-15	36-38
RHA506.3	Large	15-16	38-41
RHA506.4	X-Large	16-18	41-46
RHA506.5	XX-Large	18-20	46-51