

## Cast Shoe

Product Model : **RHA706**




**Category: Medical Equipment, Rehabilitation Supports, Fracture Aids,**

### Summary :

Looking for a comfortable and stylish way to walk while wearing a cast? Look no further than our Cast Shoes! Made from top-quality materials, these shoes are perfect for anyone who needs to wear a cast. Plus, our Walking Cast Boot makes it easy to get around town while your injury heals. Contact us today to learn more! Send your Bulk Inquiry Today.

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

**Item Code:-**

RHA706

**Features:**

The rocker bottom design assists normal gait & reduces impact coming during a heel strike. The aerospace material used on the sides makes it soft and skin friendly. The sole is specially designed to be a non-slit one. It is ambidextrous in nature.

**Applications:**

Post Cast application on foot,  
when patient is on foot cast

Cast Shoe	
Size	Item Code
Small	RHA706.1
Medium	RHA706.2
Large	RHA706.3