

## **Rib Belt**

Product Model : **RHA543**



**Category: Medical Equipment, Rehabilitation Supports, Fracture Aids,**

### **Summary :**

The Rib Belt is an orthopedic device used to provide support and compression to the rib cage. It is often prescribed following surgery or injury to the ribs or chest. The belt helps reduce pain and swelling and speeds up the healing process. The Rib Belt is made of lightweight, breathable material, and is adjustable to ensure a comfortable fit.

### **Specification :**

#### **Features:**

- Soothing soft foam laminate pads across the rib cage
- Hook and pile closure for a customizable fit
- Easy, comfortable and inconspicuous to wear
- Splint pad is placed over the site of rib fracture and the belt is wrapped around the thorax region. Fasten the hook and loop closure at a comfortable compression level.

#### **Application:**

- Rib fractures.

- Fracture of sternum.
- Cardio thoracic surgery.

Size	Size (Cms)	Size (Inches)	Item Code
S	70 to 80	27.6 to 31.5	RHA543.1
M	80 to 90	31.5 to 35.4	RHA543.2
L	90 to 100	35.4 to 39.4	RHA543.3
XL	100 to 110	39.4 to 43.3	RHA543.4
XXL	110 to 120	43.3 to 47.2	RHA543.5

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com