

## Five Functional ICU Bed Electric, Wooden Type ABS Panel

Product Model : **HF-102**



**Category: Medical Equipment, hospital furniture, ICU Beds,**

### Summary :

GWS is proud to present our Five Functional ICU Bed Electric, Wooden Type ABS Panel. This ICU bed is perfect for hospitals, emergency rooms, and other medical facilities. It features five functions: height adjustment, backrest adjustment, knee rest adjustment, Trendelenburg and reverse Trendelenburg positions, and a wooden type ABS panel. The bed is made of high-quality materials and construction that is built to last. Order your Five Functional ICU Bed Electric today!

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- ABS Head and Foot Panel with Locking and Corner Buffer.
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 460-700H mm
- Functions by Hand Remote:
  - (a) Hi-low (b)Backrest elevation (c)Kneerest elevation (d)Trendlenburg / Reverse Trendlenburg