

## Five Functional ICU Bed Mechanical (ABS Panel & Side Railing)

Product Model : **HF-112**



**Category: Medical Equipment, hospital furniture, ICU Beds,**

### Summary :

Looking for a reliable and five functional ICU bed? Look no further than the GWS Brand Five Functional ICU Bed Mechanical (ABS Panel & Side Railing) - HF112. This high-quality hospital bed is perfect for use in emergency rooms, intensive care units, and other medical settings. It features an ABS panel and side railing to provide extra security and durability. Plus, it's easy to operate and adjust, making it perfect for any setting. Order your Five Functional ICU Bed Mechanical today!

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Four Individual Side Railing, ABS Moulded.
- ABS Head and Foot Panel with Locking and Corner Buffer.
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 625-820H mm
- Functions by Smooth Crank System: (a) Hi-low (b) Backrest elevation (c) Kneerest elevation (d) Trendlenburg / Reverse Trendlenburg

### Optional:

- Four Section Mattress.