

Five Functional ICU Bed Mechanical

Product Model : **HF-113**




Category: Medical Equipment, hospital furniture, ICU Beds,


Summary :

If you're looking for a top-quality ICU bed, look no further than our GWS Brand Five Functional ICU Bed Mechanical (HF113). This hospital bed is perfect for use in emergency situations or in hospitals and clinics. It's made of durable materials and features five functionalities that make it ideal for intensive care patients. Plus, it's available at a price you'll love. Order yours today!

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com

 **Corp Office:** 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

Five Functional ICU Bed Mechanical (S.S Collapsible/Detachable Side Railing)

- Two Individual Detachable/Collapsible S.S Side Railing.
- ABS Head and Foot Panel with Locking and Corner Buffer.
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 625-820H mm

Functions by Smooth Crank System: (a) Hi-low (b) Backrest elevation (c) Kneerest elevation (d) Trendlenburg / Reverse Trendlenburg

Optional: Four section mattress