

ICU Bed Mechanical Laminated Panels

Product Model : **HF-114**




Category: Medical Equipment, hospital furniture, ICU Beds,


Summary :

GWS Brand **ICU Bed** are made with the highest quality international standards, **HF-114 ICU Bed Mechanical (Laminated Panels & S.S Side Railing)** help in recovery of a patient in the ICU and are made with much more safe and stable process, which helps motivating for patient and healthcare system.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com

 **Corp Office:** 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

ICU Bed Mechanical (Laminated Panels & S.S Side Railing)

- Two Individual Detachable/Collapsible S.S Side Railing.
- Stainless Steel Head & Foot Panel with Laminated Board & Corner Buffer
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors, 2 with Brakes.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 625-820H mm
- Functions by Smooth Crank System: (a) Hi-low (b)Backrest elevation (c)Kneerest elevation (d)Trendlenburg / Reverse Trendlenburg

Optional: Four Section Mattress.