

Examination Couch with Drawer & Cabinet

Product Model : **HF151**



Category: Medical Equipment, hospital furniture, Examination Couch & Tables,

Summary :

GWS brand examination couch is perfect for any medical facility. It features a drawer and cabinet for easy storage, and it's made from durable materials that can withstand heavy use. Plus, it's comfortable and easy to operate, making it the perfect choice for any doctor or nurse. Send your Bulk Inquiry Today!

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Made of CRCA Sheets.
- Backrest adjustable with pneumatic pump.
- Padded with foam top with high quality washable rexene.
- Provided with three drawers and three utility cabinet.
- Foldable BP apparatus tray on the side.
- Inbuilt foot step.
- Provision of saline stand.
- Pre-treated and epoxy powder coated.
- Dimension: 1870L x 560W x 810H (mm)