

## Cotton Crepe Bandage

Product Model : **SD290 & SD291**



**Category: Medical Equipment, Medical Disposables, Tapes, Bandages, Dressings & Absorbent's,**

### Summary :

A crepe bandage is a versatile medical wrap used to create a localized pressure for the treatment of muscle sprains and strains to reduce blood flow and apply an even stable pressure to restrict swelling at the site of injury. Our crepe bandages are made from a cotton-polyester latex free compound which provides adequate compression and elasticity. This improves the ease of use and prevents any reactions which may be caused due to latex. The crepe bandages come in two dimensions: 10cm x 4.5mtrs stretched length and 7.5cm x 4.5mtrs stretched length. They are widely applicable in orthopedics, sports medicine and surgical dressing.

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

Cotton crepe bandages are the perfect choice for professionals who need a bandage that is strong, yet gentle enough to not cause irritation. Our cotton crepe bandages are made from high-quality materials and are available at a competitive price. Available in Different Sizes.

<b>Size</b> (cms x mtrs.)	<b>Codes</b>
5cm x 4m	SD291.1
7.5cm x 4m	SD291.2
10cm x 4m	SD291.3
15cm x 4m	SD291.4
8cm x 4m	SD291.5
8cm x 5m	SD291.6
10cm x 5m	SD291.7