

Absorbent Pads

Product Model : **AP17**




Category: Medical Equipment, Medical Disposables, Tapes, Bandages, Dressings & Absorbent's,

Summary :

Are you looking for a quality absorbent pad? Look no further than GWS Surgicals LLP. Our disposable bed pads are made from top-quality materials designed to absorb moisture and protect your bedding from leaks and accidents. Whether recovering from surgery or just having occasional bathroom accidents, our absorbent pads will keep you comfortable and safe. Plus, our products are manufactured in India using the latest technology and strict quality controls, so you can be sure you're getting the best possible bulk quotes from underpads manufacturers in India.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

?

GWS - AbsoPads, premium bed underpads absorb up to 6 cups of liquid, so you can rest assured your sheets and mattress will remain clean and dry, night after night. Ideal alternative to incontinence diapers and other bulky incontinence aids.

?

Back sheet - White/Blue/Green PE Film

?

Top sheet - Soft non

?

woven surface

?

Super absorbent polymer (SAP) ensures complete dryness and comfort.

?

Diamond embossed pattern for uniform and faster distribution of liquids.

?

Also available with silicone release tape for back sheet sticking.

?

Available in Packing: 1pc. | 3pcs. | 5pcs. | 10pcs.

Size	Size (in cms)	Item Code
Small	60cms x 40cms	AP17.1
Medium	60cms x 60cms	AP17.2
Large	60cms x 90cms	AP17.3

X-Large	60cms x 180cms	AP17.4
---------	----------------	--------

Product MOQ:

50,000 pcs.

Master Carton (Dimension):-

60cms x 42cms x 40cms

Qty Inside Each Master Carton: -

200pcs. / 20packs of 10pcs.