

DuraLOCK Small Fragment Locking Implant Set

Product Model : **SGL14**



Category: Orthopaedic Implants, Orthopaedic Sets,

Summary :

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

SGL14 duraLOCK Small Fragment Implant Set, S.S.		
	Set Consists of : (Plates & Screws)	Qty
SGL14.1	DuraLOCK Plate 3.5 mm - S.S.	
	5, 6, 7, 8, 9, 10, 11, 12, 14 Holes (2 each)	16
SGL14.2	DuraLOCK T-Plate with 3 Head Holes 3.5 mm - S.S.	
	3, 4, 5, 6, 7 Holes (2 each)	10
SGL14.3	DuraLOCK T-Plate with 4 Head Holes 3.5 mm - S.S.	
	3, 4, 5, 6, 7 Holes (2 each)	10
SGL14.4	DuraLOCK One Third Tubular Plate 3.5 mm - S.S.	
	5, 6, 7, 8, 9, 10, 12, 14 Holes (2 each)	16
SGL14.5	DuraLOCK Proximal Humerus Plate 3.5 mm - S.S.	
	3, 5, 7 Holes (2 each)	6
SGL14.6	DuraLOCK T-Plate Oblique 3.5 mm - S.S. - Left/Right	
	3, 4, 5, 7 Holes (2 each) - Left	8
	3, 4, 5, 7 Holes (2 each) - Right	8
SGL14.7	DuraLOCK Reconstruction Plate 3.5 mm - S.S.	
	5, 6, 7, 8, 10, 12, 14 Holes (2 each)	14
SGL14.8	DuraLOCK Self Tapping Screw 2.7 mm, Hexagonal Socket - S.S.	
	Length: 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm (3 each)	42
SGL14.9	Cortex Screw 2.7 mm, Hexagonal Socket - S.S.	
	Length: 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm (3 each)	42
SGL14.10	DuraLOCK Self Tapping Screw 3.5 mm- S.S.	
	Length: 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 45mm, 50mm, 55mm, 60mm (4 each)	80
SGL14.11	Cortex Screw 3.5 mm, Hexagonal Socket - S.S.	
	Length: 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 45mm, 50mm, 55mm, 60mm (4 each)	80
SGL14.12	Cancellous Screw 4.0 mm, Fully Threaded, Hexagonal Socket - S.S.	
	Length: 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 45mm, 50mm, 55mm, 60mm (3 each)	60