

Cannulated Headless Compression Screw Instrument Set 3.0mm/3.5mm/4.0mm/5.0mm

Product Model : **HS50**



Category: Orthopaedic Implants, Orthopaedic Sets,

Legacy of Excellence in Healthcare...

Specification

HS50	Cannulated Headless Compression Screw Instrument Set 3.0mm/3.5mm/4.0mm/5.0mm	Qty
	Set Consists of:	
1	Clean Wire Ø1.5*200	1
2	Canuated Drill Bit Ø3.5*Ø1.5*200 with Limitator	1
2.1	Canuated Drill Bit Ø3.5*Ø1.5*200	1
2.2	Ø3.5 Limitator	1
3	Guider Wire 1.5*200	1
4	Guider Wire 1.5*200	1
5	Cannulated Countersink Ø5.0*Ø1.5*150	1
6	Hex Screwdriver SW3.0	1
7	Cannulated Hex Screwdriver SW3.0	1
8	Wire Sleeve Ø1.6	1
9	Drill Sleeve Ø3.5	1
10	Canuated Drill Bit Ø2.8*Ø1.2*150 with Limitator	1
10.1	Canuated Drill Bit Ø2.8*Ø1.2*150	1
10.2	Ø2.8 Limitator	1
11	Guider Wire 1.2*150	1
12	Guider Wire 1.2*150	1
13	Cannulated Countersink Ø3.5*Ø1.3*150	1
14	Hex Screwdriver SW2.5	1
15	Cannulated Hex Screwdriver SW2.5(Ø1.3)	1
16	Wire Sleeve Ø1.3*3.0	1
17	Drill Sleeve Ø2.8	1
18	Drill Sleeve Ø2.0	1
19	Wire Sleeve Ø1.1*2.0	1
20	Canuated Drill Bit Ø2.0*Ø1.1*150 with Limitator	1
20.1	Canuated Drill Bit Ø2.0*Ø1.1*150	1
20.2	Ø2.0 Limitator	1
21	Cannulated Hex Screwdriver SW2.0	1
22	Hex Screwdriver SW2.0	1
23	Guider Wire 1.0*150	1
24	Guider Wire 1.0*150	1
25	Guider Wire 1.0*150	1
26	Guider Wire 1.0*150	1
27	Cannulated Countersink Ø3.0*Ø1.1*150	1
28	Drill Sleeve Ø2.5	1
29	Wire Sleeve Ø1.1*2.5	1
30	Canuated Drill Bit Ø2.5*Ø1.1*150 with Limitator	1
30.1	Canuated Drill Bit Ø2.5*Ø1.1*150	1
30.2	Ø2.5 Limitator	1
31	Clean Wire Ø1.2	1
32	Depth Gauge 10-60mm&10-70mm	1
33	Hex Key SW2.5	1
34	Screw Holding Forcep	1
35	Quick Coupling Handle	1