

Spoon Plate, 4.5mm

Placa Cuchara 4.5mm

Model: 1315

Categories: Orthopaedic Implants, Large Fragment

Spoon Plate, 4.5mm
Page not found or type unknown

Summary

The **Spoon Plate 4.5mm** is an orthopaedic implant available in Stainless Steel 316L and Titanium, designed for fracture fixation within the 4.5mm Large Fragment System. It features a contoured spoon-shaped head, multiple countersunk screw holes, and an oblong slot for positioning and adjustment during fixation.

Specifications

Material: Stainless Steel 316L & Titanium

System: 4.5mm Large Fragment

Design: Spoon-shaped plate with widened head

Surface Finish: Smooth, polished

Hole Configuration: Multiple countersunk screw holes and oblong adjustment hole

Description

The **Spoon Plate 4.5mm**, manufactured by **GWS Surgicals LLP**, is an orthopaedic implant available in **Stainless Steel 316L and Titanium**, designed for use in fracture fixation procedures. The plate features a widened, spoon-shaped head with a contoured profile to support stable bone fixation. It incorporates multiple countersunk screw holes for secure screw placement and an oblong hole for positioning and adjustment during surgery. Manufactured from high-quality materials, the plate offers strength, durability, and corrosion resistance. Its smooth, polished surface and precision-engineered design make it a suitable component of the **4.5mm Large Fragment Orthopaedic Implant System**.