

# Hospital Instrument Trolley with Two Drawers

Carro Hospitalario para Instrumental con Dos Cajones y Bandejas de Acero Inoxidable

**Model: HF231**

Categories: Hospital & Medical Furniture, Hospital Trolleys

Hospital Instrument Trolley with Two Drawers

## Summary

Hospital Instrument Trolley with Two Drawers (HF231) features a stainless steel top tray, two storage drawers, a lower shelf, and a side accessory holder. Designed for convenient storage and movement of medical instruments and supplies in hospitals, clinics, and healthcare facilities.

## Specifications

| Parameter      | Details                                      |
|----------------|----------------------------------------------|
| Product Name   | Hospital Instrument Trolley with Two Drawers |
| Item Code      | HF231                                        |
| Product Type   | Medical Instrument Trolley                   |
| Top Tray       | Stainless Steel                              |
| Drawers        | 2                                            |
| Drawer Handles | Metal                                        |
| Lower Shelf    | Stainless Steel                              |

|                         |                                                       |
|-------------------------|-------------------------------------------------------|
| Frame                   | White-finished metal                                  |
| Side Accessory Holder   | Perforated holder                                     |
| Guard Rails             | Stainless steel tubular rails                         |
| Mobility                | Swivel castor wheels                                  |
| Wheel Brakes            | Visible on some castors                               |
| Application             | Medical instrument and supply storage                 |
| Usage Areas             | Hospitals, clinics, nursing homes and treatment rooms |
| Brand                   | GWS                                                   |
| Manufacturer / Supplier | GWS Surgicals LLP                                     |

## Description

The Hospital Instrument Trolley with Two Drawers (HF231) is designed for convenient storage, organization, and transportation of medical instruments and supplies.

It features a stainless steel top tray with protective guard rails, two spacious drawers with metal handles, and a lower stainless steel shelf for additional storage. A perforated side holder provides convenient space for frequently used accessories.

The trolley is supported by a sturdy, white-finished frame and fitted with swivel castor wheels for easy movement within healthcare facilities.

Suitable for hospitals, nursing homes, clinics, treatment rooms, and medical departments.