
5+3 Twin

Model: OTE53

Categories: Medical Equipment, Operation Light, Elite Series Operation Light



Description

Technical Data	
Intensity	1,60,000 + 1,00,000 Lux ± 10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Reduction Index	93 RA
LED Life	50000 Hrs.
Number of LED	105 + 63
Diameter of Light	600mm x 2
Brightness Control	Capacitive Touch Panel

Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional