Foley Balloon Catheter Siliconized

Sonda Foley con Balón Siliconizada UD761

Model: UD761

Categories: Medical Equipment, Medical Disposables, Urology



Summary

The GWS Foley Balloon Catheter UD761 is a siliconized latex urinary catheter with atraumatic drainage eyelets and a symmetrical retention balloon, available in sizes 8–26 FR. It is single-use, ethylene oxide sterilized, and designed for continuous bladder drainage in healthcare settings.

Specifications

Material: Siliconized Latex

Surface: Smooth, Siliconized

Drainage: Atraumatic Drainage Eyelets

Balloon: Symmetrical Retention Balloon

Inflation Port: Colour-Coded Valve

Connector: Universal Funnel Connector

Available Sizes: 8 FR to 26 FR

Sterilization: Ethylene Oxide (EO)

Usage: Single Use Only

Packaging: 10 Units per Box

Manufacturer: GWS Surgicals LLP

Country of Origin: India

Description

The GWS Foley Balloon Catheter UD761 is a balloon retention urinary catheter manufactured from siliconized latex for bladder catheterization and continuous urine drainage. The smooth siliconized surface, atraumatic drainage eyelets and symmetrical retention balloon support safe catheter positioning and drainage. A colour-coded inflation valve enables size identification, while the universal funnel connector accommodates compatible urine collection systems. Available in ten French sizes from 8 FR to 26 FR, the catheter is ethylene oxide sterilized and supplied as a single-use disposable device. It is designed for use by qualified healthcare professionals in hospitals, clinics and other appropriate healthcare settings.