Dual Power Suction Unit (Electric cum Manual)

Model: SU-21

Categories: Medical Equipment, Suction Machine, Electric cum Manual Suction Units



Description

- Housing: MS Powder Coated, SS Top on two wheels
- Capacity (Electrical): -600 mm Hg ± 10 , 45 Ltrs/ min
- Capacity (Pedal): -600 mm Hg ± 10 , 200ml/ Push Stroke
- Pump Type: Electric – Diaphragm / Pedal – Piston
- Jars: 2 x 2 Ltrs Polycarbonate Jar
- Filter: Bacterial Filter Autoclavable/Reusable
- Tubing: 10mm ID x 2 mtr. (PVC)
- Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration.
- Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request)
- Noise Level: 55 dB A ± 3
- Dimension & Weight: 42 x 39 x 64 cms, 15.5 Kg
- Overflow Safety: Dual Mechanical